

# Model Z Material Best Practice



**Desktop Health™**

**Manufacturer:**  
EnvisionTEC US LLC  
A Desktop Metal Company  
15162 S. Commerce Dr. Dearborn,  
Michigan 48120  
USA

# Table of Contents

<b>About Model Z</b> .....	<b>3</b>
Identification .....	3
Applicable Printers.....	3
<b>Getting Started</b> .....	<b>4</b>
Primary Supplies .....	4
Capture Patient Data .....	4
Design Models for Model Z.....	5
<b>Software</b> .....	<b>6</b>
Orient Models Envision One RP Software .....	6
Autopilot   Models .....	6
<b>Print Preparation</b> .....	<b>7</b>
Mix Material .....	7
Fill Material Tray .....	7
Print with Model Z Material.....	7
<b>Post-Processing</b> .....	<b>8</b>
Clean Printed Models .....	8
Dry Models .....	9
Post Cure Printed Models .....	9

Model Z material was developed for use in the dental industry as dental models. The manufacturer does not accept any liability for damages caused by any other use of the material. Furthermore, before using the material, the user must independently check for its suitability and applicability for the intended use, especially if not specified in this document. Einstein is a registered trademark of Desktop Metal, Inc. Envision One is a registered trademark of EnvisionTEC GmbH. Desktop Health is a registered trademark of Desktop Metal, Inc. Vaseline is a registered trademark of Conopco, Inc. Siliform is a registered trademark of Keilhauer Keilhauer Industries II Ltd.

# About Model Z

## Identification

Model Z is as effective as it is fast. This low viscosity model material is ideal for thermoforming models, diagnostic models and orthodontic appliance models requiring a smooth surface resolution, easy clean-up and a superior fit and finish.


This technical guide details the best practices for preparing models, post-processing, and finishing.

## Applicable Printers

This material is tested and approved for the following printers:

- Einstein™
- Envision One cDLM
- D4K Pro

Fig. 1 EINSTEIN™ FRONT VIEW



Model Z material was developed for use in the dental industry as dental models. The manufacturer does not accept any liability for damages caused by any other use of the material. Furthermore, before using the material, the user must independently check for its suitability and applicability for the intended use, especially if not specified in this document. Einstein is a registered trademark of Desktop Metal, Inc. Envision One is a registered trademark of EnvisionTEC GmbH. Desktop Health is a registered trademark of Desktop Metal, Inc. Vaseline is a registered trademark of Conopco, Inc. Siliform is a registered trademark of Keilhauer Keilhauer Industries II Ltd.

# Getting Started

## Primary Supplies

The following supplies are required to print Model Z material:

- 99% isopropyl alcohol (IPA).
- Air compressor.
- Cone-shaped paint filter, *Starter Kit item*.
- Curing unit: Otofash SAP Part # ACC-00-0007, or PCA 4000 SAP Part # ACC-06-1000.
- Dual Motion Bottle Roller, SAP Part # ACC-26-1000 (110V) and ACC-26-1000 (220V).
- Nitrile gloves.
- Paint scraper, *Starter Kit item*.
- Paper towels.
- Plastic funnel.
- Rubber spatula, *Starter Kit item*.
- Spray bottle with 99% IPA.
- Snips, precision blade, or similar tool.
- Storage container for material, sealable and opaque.
- Washing unit: PWA 2000, SAP # ACC-22-2000.

## Capture Patient Data

A digital impression can be accomplished with a handheld intraoral scanner and CBCT scan, or with a traditional impression and a desktop box scanner.


Envision One RP Software is compatible with the universal .STL file format and is thus compatible with almost all dental CAD and model design software as well as digital design services. Models may be designed in-house or outsourced to a design partner.

Model Z material was developed for use in the dental industry as dental models. The manufacturer does not accept any liability for damages caused by any other use of the material. Furthermore, before using the material, the user must independently check for its suitability and applicability for the intended use, especially if not specified in this document. Einstein is a registered trademark of Desktop Metal, Inc. Envision One is a registered trademark of EnvisionTEC GmbH. Desktop Health is a registered trademark of Desktop Metal, Inc. Vaseline is a registered trademark of Conopco, Inc. Silifform is a registered trademark of Keilhauer Keilhauer Industries II Ltd.

## Design Models for Model Z

Hollow dental models printed in Model Z must have a minimum wall thickness of 3 mm, *Fig. 2*.

**Fig. 2** HOLLOW MODEL WALL THICKNESS



It is recommended to add channels or drainage holes to hollow models. This allows uncured material to drain from the hollow feature during the printing process.

# Software

## Orient Models Envision One RP Software


Orient models in Envision One RP Software with the flat base side down, facing the build platform.

- **Spacing:** place models a minimum of 2.5 mm apart.
- **Level at build platform:** place models 0 mm from the build platform.
- **Resolution:** 100 µm Z resolution.

## Autopilot | Models

**Autopilot | Models** in Envision One RP Software can be used to quickly orient all selected orthodontic models.

Fig. 3 AUTOPILOT | MODELS



**Models:** Aligns the largest flat surface with the build platform, spaces the models across the build platform, and transfers the job automatically to the printer.

# Print Preparation

## Mix Material

Model Z material must be mixed in the material bottle prior to use:

1. Place the sealed material bottle on the Dual Motion Bottle Roller for a minimum of 30 minutes.
2. Wait for bubbles to subside before filling the material tray.
3. Mix material in the material tray gently with the rubber spatula from the Starter Kit before each print. The material should be a uniform color.

Ensure there are no small cured particles in the material. If found, then the material must be filtered using the plastic funnel, cone-shaped paint filter, and a spare material bottle. *See the [Knowledge Base](#) for filtering instructions.*

## Fill Material Tray

Do not overfill the material tray. Overfilling can cause the material to overflow when the build platform moves down at the start of the print job.

To add more material to the printer, carefully pour material into the material tray between prints. Adding material while the print is paused, or during a print, will cause a small shift line in the model. *See the [Knowledge Base](#) for instructions adding material.*

## Print with Model Z Material

To start the print, follow instructions in the printer's User Manual.

To remove the models from the build platform after the print is complete, follow instructions in the printer's User Manual. *See the [Knowledge Base](#) for the latest User Manual.*

# Post-Processing


## Clean Printed Models

The PWA 2000 is the recommended parts washer, *Fig. 4*. Always wear gloves when handling uncured material and alcohol.



**Important: Do not expose Model Z to alcohol for longer than 5 minutes.** Excess exposure to alcohol will dry out the model.

**Fig. 4** PWA 2000 FRONT VIEW



### Clean models the PWA 2000:

1. Open the washing compartment lid.
2. Lift the handle to raise the interior grate to the highest position.
3. Fill the washing compartment with 99% IPA.
4. Place the model on the grate and gently lower the handle to submerge the model in 99% IPA.
5. Close the washing compartment lid and lock in place.
6. Plug in the power cable to turn on the PWA-2000.
7. Using the touchscreen, select the **High** washing program. Set the timer to 00:03:00, or 3 minutes. Press **Start**.  
→ *The PWA 2000 will immediately begin the set washing cycle.*
8. Remove the model as soon as the program is complete.
9. Spray the models with the spray bottle filled with 99% IPA.
10. Use compressed air to remove all IPA from the surface of the model as soon as possible.

Model Z material was developed for use in the dental industry as dental models. The manufacturer does not accept any liability for damages caused by any other use of the material. Furthermore, before using the material, the user must independently check for its suitability and applicability for the intended use, especially if not specified in this document. Einstein is a registered trademark of Desktop Metal, Inc. Envision One is a registered trademark of EnvisionTEC GmbH. Desktop Health is a registered trademark of Desktop Metal, Inc. Vaseline is a registered trademark of Conopco, Inc. Silifform is a registered trademark of Keilhauer Keilhauer Industries II Ltd.



## Dry Models

Models must be completely dry before post curing –

1. Place the models on a clean paper towel lined surface.
2. Air dry in ambient room temperature / humidity for 10 min.

## Post Cure Printed Models



Cure the models using one of the following methods:

**Otoflash:** 2 cycles for 500 flashes, flip models between cycles.

See the [Knowledge Base](#) for instructions setting an Otoflash curing cycle.

**PCA 4000:** 1 Minute - 20° C - 100% Power.

See the [Knowledge Base](#) for instructions setting a PCA 4000 curing cycle.



Place models into the curing unit with as much space between models as possible. Models should never touch one another while curing. Let models cool completely before handling them or starting the next cycle. Flip models between cycles for an even cure.



**Important:** Desktop Health does not support third-party curing units.

## Release Agents

Desktop Health does not recommend a specific isolation or release agent for the production of orthodontic appliances. Vaseline®, Siliform®, or any other release agent may be used. Follow instructions from the manufacturer.

Model Z material was developed for use in the dental industry as dental models. The manufacturer does not accept any liability for damages caused by any other use of the material. Furthermore, before using the material, the user must independently check for its suitability and applicability for the intended use, especially if not specified in this document. Einstein is a registered trademark of Desktop Metal, Inc. Envision One is a registered trademark of EnvisionTEC GmbH. Desktop Health is a registered trademark of Desktop Metal, Inc. Vaseline is a registered trademark of Conopco, Inc. Siliform is a registered trademark of Keilhauer Keilhauer Industries II Ltd.